

MODERN SURGICAL COMPANY

The world Of Ophthalmic Equipments
ADDRESS: 101/C, C.R. AVENUE KOLKATA- 700073
+91-9830302017/9308241770



SPECIFICATION

SL.NO	Parameter	<i>EOSTER SA- SEMI AUTOMATIC</i>
1	Chamber Size	35 ltrs
2	Chamber Outer Body	Stainless- SS304 MS
3	Sterilization Agent	Standard -EO gas cartridges (100%)
4	Sterilization Cycle	Warm Cycle - 55°C
5	Control Process	Automated, PLC
6	Aeration Method	Aeration is done inside chamber by means of vacuum pump, Inbuilt Aerator.
7	Vacuum Pump	Dry Vacuum Pump
8	Heater	Industrial heater uniformly heats the chamber, 1400 W
9	Power	1 ph, 230V stabilized Power
10	External Dimension	28"X28"X18" Approx
11	Weight	60 Kgs approx
12	Sealing machine	One
13	ETO Gas	25 pcs
14	Packing Paper	2 rolls
15	Trolley Stand	2 pc stainless
16	ETO Indicator paper	1 pkt